

Colorectal Air Syringes

“Routine testing of colorectal anastomoses reduces the postoperative leak rate.”
*American Society of Colon and Rectal Surgeons**

**STANDARDIZE THE
INTRAOPERATIVE
AIR LEAK TEST**



COLORECTAL AIR SYRINGE-LP

Inflates to 30–40 cmH₂O above intraperitoneal pressure when the abdomen is insufflated to 15 mmHg

White for LAPAROSCOPIC and ROBOTIC surgery
Catalog No. ALT-50

- Occlude the anus
- Regulate inflation pressure
- Signal at targeted pressure
- Sterile, compact, and disposable
- Rectal tube for bowel deflation included



COLORECTAL AIR SYRINGE-OP

Inflates to 30–40 cmH₂O

Black for OPEN surgery
Catalog No. ALT-30

U.S. Patent Nos. 8,216,159 and 10,123,712

*Feingold D, Steele SR, Lee S, Kaiser A, Boushey R, Buie WD, Rafferty JF. Clinical Practice Guideline Task Force of the American Society of Colon and Rectal Surgeons. Practice parameters for the treatment of sigmoid diverticulitis. Dis Colon Rectum 2014;57:284-294

TfS Tools
for
Surgery

www.toolsforsurgery.com/ALT

Tools for Surgery, LLC
8 Technology Drive, #100 • East Setauket, NY 11733
P: 631-444-4448 • F: 631-689-2943
info@toolsforsurgery.com